

---

**Monoclonal Anti- Brain natriuretic peptide, BNP (Detection Ab)****Catalog#** BDA1029**Lot #** Check on the product label**Size** 1 mg(5.28 mg/ml)**Isotype** IgG2b**Clone #** 2D6**Host** Mouse**Reactivity**

Human

**Product Form** Liquid**Purification & Buffer**

Protein A or G purified and supplied in 0.01 M PBS (pH7.4) without preservative.

**Purity** >95% by HPLC & SDS-PAGE**Immunogen** Recombinant BNP protein.**Recommend Application**

ELISA

Chemiluminescent immunoassay, CLIA

Flow Fluorescence Immunmicrobeads

Assay, FFIA

Other applications have not been tested.

The optimal dilutions should be determined by end user.

**Matched antibody pair**

Capture Ab: BNP mAb (clone # B4)

Detection Ab: BNP mAb (clone # 2D6)

**Storage Instruction**

Aliquot and store at -20°C for long term (at least one year).

Avoid repeated freeze and thaw cycles.

**Background**

Brain natriuretic peptide or B-type natriuretic peptide (BNP) (also ventricular natriuretic peptide or natriuretic peptide B) is a 32-amino acid polypeptide secreted

by the ventricles of the heart in response to excessive stretching of heart muscle cells (cardiomyocytes). BNP is secreted attached to a 76 - amino acid N-terminal fragment in the prohormone called NT-proBNP (BNPT), which is biologically inactive. BNP is synthesized as a 134-amino acid preprohormone (preproBNP), encoded by the human gene NPPB.

**Reference**

1. Ziskoven D, Forssmann WG, Holthausen U, Menz G, Addicks K, Rippegater G (1989). "Calcium Calmodulin antagonists Influences the release of Cardiodilatin/ANP from Atrial Cardiocytes". In Kaufmann W, Wambach G. Handbook Endocrinology of the Heart. Berlin: Verlag: Springer. pp. 233–4.

---

**FOR RESEARCH USE ONLY, NOT FOR DIAGNOSTIC AND CLINICAL USE.**

Chongqing Biospes Co., Ltd Tel: +86-23-67567091 Fax: +86-23-67745923

7F, Bldg B, High-tech Venture Park, # 107 Erlang Chuangye Rd, Jiulongpo District, Chongqing, 400039, China

[www.biospes.com](http://www.biospes.com)