
SDS-PAGE Electrophoresis Buffer

Catalog# BWR1030

Size: 1 L

Lot # Check on the product label

Storage Instruction: Store at 2-8°C for 6 months.

Content: 100mM Tris, 100mM tricine, 0.1% SDS.

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC AND CLINICAL USE.

Chongqing Biospes Co., Ltd Tel: +86-23-67567091 Fax: +86-23-67745923
7F, Bldg B, High-tech Venture Park, # 107 Erlang Chuangye Rd, Jiulongpo District, Chongqing, 400039, China
www.biospes.com