
Coomassie Brilliant Blue G250

Catalog# BWR1043

Size: 5g

Lot # Check on the product label

Storage: Store at room temperature.

Formula: C₄₇ H₄₈ N₃ Na O₇ S₂ (sodium salt)

CAS number: 6104-58-1

Molecular weight: 854.02

Solubility: Slightly soluble in water, best to dissolve first in methanol or ethanol.

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC AND CLINICAL USE.

Chongqing Biospes Co., Ltd Tel: +86-23-67567091 Fax: +86-23-67745923

7F, Bldg B, High-tech Venture Park, # 107 Erlang Chuangye Rd, Jiulongpo District, Chongqing, 400039, China

www.biospes.com