

---

**Coomassie Brilliant Blue R250**

**Catalog#** BWR1044

**Size:** 5 g

**Lot #** Check on the product label

**Storage:** Store at room temperature.

**Formula:**  $C_{45}H_{44}N_3NaO_7S_2$  (sodium salt)

**CAS number:** 6104-59-2

**Molecular weight:** 825.97

**Solubility:** Essentially insoluble in water, dissolve first in methanol.

---

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC AND CLINICAL USE.

Chongqing Biospes Co., Ltd Tel: +86-23-67567091 Fax: +86-23-67745923  
7F, Bldg B, High-tech Venture Park, # 107 Erlang Chuangye Rd, Jiulongpo District, Chongqing, 400039, China  
[www.biospes.com](http://www.biospes.com)